



Philips Respironics V200 Esprit Oxygen Inlet Filter Replacement

Use this guide to replace the oxygen inlet filter in the Philips Respironics V200 Esprit.

Written By: Grace Gius




Figure 8-2: Front Panel — Text Version

INTRODUCTION

This guide is for the Philips Respironics V200 / Esprit Ventilator, Part Number 580-1000-02 H. The oxygen water trap/inlet filter assembly consists of a 5-micron (μ) filter to remove particulate (both dry and liquid) from the oxygen gas supply, a bowl with drain for accumulated water, and an oxygen inlet connector (Chapter 3.1.9 on page 3-4 of [document|7525|Service Manual]).

Information in this guide was sourced from Chapter 10.1.3 on page 10-5 of the Respironics V200/Esprit Ventilator [document|7525|Service Manual]. Further information on general filter replacements for this Ventilator can be found in Chapter 9.1 on pages 9-3 and 9-4 of the [document|7525|Service Manual].



PARTS:

- [Filter](#) (1)
5-micron
-

Step 1 — Oxygen Inlet Filter




Figure 8-2: Front Panel — Text Version

- Turn off the device.
- On/Off switch located on the Front Panel as seen in Figure 8-2 (from the [Operator Manual](#)).

Step 2



- Disconnect from AC power.

Step 3

- Disconnect the oxygen hose.
- Oxygen hose is connected to the Inlet Connector (section 14 of Figure 5-1, diagram from [Operator Manual](#)).




Figure 5-1: Back Panel

Step 4





Figure 5-1: Back Panel

- Unscrew the water trap bowl (page 10-5 in [Service Manual](#)).
- Located on the back panel of the ventilator (section 15 of Figure 5-1 from [Operator Manual](#)).

Step 5



- Separate the filter element from the black plastic retainer (page 10-5 of [Service Manual](#)).

Step 6



- Attach new filter element to the black plastic retainer (page 10-5 of [Service Manual](#)).

Step 7



- Screw on water trap bowl (page 10-5 of [Service Manual](#)).

Step 8



- Reconnect the oxygen hose.
 - According to page 2-2 of the [Operators Manual](#), care in the routing of the oxygen inlet hose should be exercised to ensure it is not exposed to mechanisms that could cause damage by cutting or heating/melting.
-
- Reconnect to AC Power.
 - To reassemble your device, follow these instructions in reverse order.